

Zefon Gas Sampling Lines

- Economical
- Disposable, Single Patient Use
- Non-Sterile
- Sizes Available for Any Monitor
- 100% Flow Tested
- Top Quality Medical Grade Tubing
- Filters & Nafion® Dryer Tube Options
- Custom Configurations Available



LATEX FREE

Ordering Information

8' PVC Gas Sampling Lines

Product #	Tube ID	Connectors
ZML-9300	.040"	Male/Male
ZML-9400	.040"	Male/Female
ZML-9330	.050"	Male/Male
ZML-9430	.050"	Male/Female
ZML-9320	.060"	Male/Male
ZML-9420	.060"	Male/Female

10' PVC Gas Sampling Lines

Product #	Tube ID	Connectors
ZML-9500	.040"	Male/Male
ZML-9600	.040"	Male/Female
ZML-9530	.050"	Male/Male
ZML-9630	.050"	Male/Female
ZML-9520	.060"	Male/Male
ZML-9620	.060"	Male/Female

8' Co-Extruded PE/PVC Gas Sampling Lines

Product #	Tube ID	Connectors
ZML-9350	.040"	Male/Male
ZML-9450	.040"	Male/Female
ZML-9340	.047"	Male/Male
ZML-9440	.047"	Male/Female

10' Co-Extruded PE/PVC Gas Sampling Lines

Product #	Tube ID	Connectors
ZML-9550	.040"	Male/Male
ZML-9650	.040"	Male/Female
ZML-9540	.047"	Male/Male
ZML-9640	.047"	Male/Female

All Zefon Gas Sampling Lines come individually packaged with 50 units per case. Call for custom configurations.

Attachment Options

Attachments are available pre-installed on any line. To order, select the option you want and adjust the product number to reflect the option. To do this, simply replace the last digit of the product number with the option number. Example: To add option #3 to a ZML-9520, change the number to ZML-9523.

Option #1 - Add a 5.0µm Filter to the sampling line

Option #3 - Add a 12" Nafion® dryer segment to the sampling line

Option #6 - Add both a 5.0µm filter and a 12" Nafion® dryer segment to the line.

Nafion is a registered trademark of Dupont

Distributed By:

