

NeoCheck™



Disposable Blood Pressure Cuffs

For over twenty years Zefon International has been designing and manufacturing superior blood pressure cuffs for leading blood pressure measurement instrument manufacturers. With a complete range of sizes, materials and connectors, Zefon brand cuffs represent the best value to bring consistent performance, patient comfort and economy to your medical setting.

Features

- Disposable; Eliminates costly disinfection procedures and contamination concerns.
- Cost Effective; Economical enough for single patient use.
- Increase infection control with single patient use cuffs.
- Available in soft white fabric material or flexible white vinyl material.
- Meets AAMI and AHA Cuff Sizing Recommendations.
- Color coded printing on each cuff makes it easy to identify cuff size.
- Latex Free.
- Designed as a direct replacement for OEM cuffs provided by monitor manufacturers.
- Full range of connectors and adapters allow standardization of cuffs across the whole hospital.
- Easy connection to electronic NIBP monitors or manual blood pressure measurement systems.




Zefon International Inc.
5350 SW 1st Lane, Ocala, FL 34474 USA
Tel: +1 3528548080; Fax: +1 3528547480
www.zefon.com

MM00003 Rev 5 2010-07

EC REP Blue Box Medical Limited
Unit 29, New Forest Enterprise Centre,
Rushington Business Park; Chapel Lane
Totton, Southampton, SO40 9LA United Kingdom
Phone: +023 8066 9000; Fax: +023 8066 9002

Ordering Information


Single Tube Cuffs



			Neonatal Monitor versions from: Philips/HP Siemens Spacelabs Datascope	Adult Monitors From: CAS 740/750 Criticare 508 Datascope Accutor Plus & Passport 2 Philips/HP Siemens Spacelabs	Adult Monitors From: BCI Colin Datascope Passport XG & Expert Nihon Kodan Lifescope 9 & 12 Welch Allyn Vital Signs/Atlas (Before July 2000)	Adult Monitors From: Welch Allyn Vital Signs/Atlas (after July 2000)	Adult Monitors From: Welch Allyn Protocol Propaq
Cuff Size	Color Code	Limb Circumference	1-Tube Male Slip Luer	1-Tube Bayonet Connector	1-Tube Female Luer Connector	1-Tube Screw Type Connector	1-Tube Female Twist Lock Connector
Neonatal #1	Orange	3-6cm	ZNC4501L	ZNC4501HP	ZNC4501F	ZNC4501D	ZNC4501MQF
Neonatal #2	Light Blue	4-8cm	ZNC4601L	ZNC4601HP	ZNC4601F	ZNC4601D	ZNC4601MQF
Neonatal #3	Green	6-11cm	ZNC4701L	ZNC4701HP	ZNC4701F	ZNC4701D	ZNC4701MQF
Neonatal #4	Dark Blue	7-13cm	ZNC4801L	ZNC4801HP	ZNC4801F	ZNC4801D	ZNC4801MQF
Neonatal #5	Red	8-15cm	ZNC4901L	ZNC4901HP	ZNC4901F	ZNC4901D	ZNC4901MQF

Product numbers listed are for NeoCheck Soft cuffs. To order vinyl NeoCheck cuffs, change the ZNC prefix to ZNCV. All cuffs are packaged with 10 cuffs per box

Double Tube Cuffs



			Neonatal Monitor versions from: Critikon Dinamap Marquette MDE	Adult Monitors From: Critikon® Dinamap® Criticare 507 Invivo RESEARCH MDE Nihon Kodan Lifescope L Spacelabs	Adult Monitors From: Datex-Ohmeda	Adult Monitors From: GE Medical Eagle 4000 Tram 100/200/300	Manual Inflation
Cuff Size	Color Code	Limb Circumference	2-Tube Male Slip Luer	2-Tube Screw Type Connector	2-Tube Male/Female Twist Lock Connector	2-Tube Female Twist Lock Connector	2-Tube Manual Inflation with Male Luer
Neonatal #1	Orange	3-6cm	ZNC4502L	ZNC4502D	ZNC4502MQMF	ZNC4502MQF	ZNC4512L
Neonatal #2	Light Blue	4-8cm	ZNC4602L	ZNC4602D	ZNC4602MQMF	ZNC4602MQF	ZNC4612L
Neonatal #3	Green	6-11cm	ZNC4702L	ZNC4702D	ZNC4702MQMF	ZNC4702MQF	ZNC4712L
Neonatal #4	Dark Blue	7-13cm	ZNC4802L	ZNC4802D	ZNC4802MQMF	ZNC4802MQF	ZNC4812L
Neonatal #5	Red	8-15cm	ZNC4902L	ZNC4902D	ZNC4902MQMF	ZNC4902MQF	ZNC4912L

Product numbers listed are for NeoCheck Soft cuffs. To order vinyl NeoCheck cuffs, change the ZNC prefix to ZNCV. All cuffs are packaged with 10 cuffs per box

Distributed By:



5350 SW 1st Lane,
Ocala, FL 34474 USA
Tel: +1 3528548080
www.zefon.com