

1 Patient 1 Cuff

Infection Control

from Admission to Discharge



The PoppaTM

Universal NIBP Cuff Connection System

A true universal connection system for NIBP Cuffs that allows standardization of connections on NIBP monitors and blood pressure cuffs.

Features:

- **Unique design** makes connecting and disconnecting easy, yet still provides a leak free seal.
- **No more splitters, Y-pieces, or other adapting type systems.**
- **Simple, easy to clean and disinfect design** with NO moving parts.
- **Improve safety** - The POPPA has purposely been designed to be completely incompatible with Luer connections to avoid any possible misconnection with intravenous lines and catheters.
- **Follows the Joint Commission recommendations*** to avoid purchasing non-intravenous devices that can physically mate with a female Luer IV line adapter.
- **Inventory management** - Stock only one style cuff for each size range that can connect to any NIBP monitor.
- **Cut costs** - A standardized connector lets you use just one disposable cuff per patient from admission to discharge.
- **Reduce hospital acquired infections** - implement a 1-patient 1-cuff policy were no one shares a cuff.

*Tubing misconnections—a persistent and potentially deadly occurrence - Joint Commission Sentinel Event Alert, Issue 36, April 3rd, 2006. www.jointcommission.org



5350 SW 1st Lane
Ocala, FL 34474 USA
Tel: +1 3528548080
Fax: +1 3528547480
www.zefon.com

EC REP


Blue Box Medical Limited
Unit 29, New Forest Enterprise Centre,
Rushington Business Park; Chapel Lane
Totton, Southampton, SO40 9LA United Kingdom
Phone: +023 8066 9000; Fax: +023 8066 9002

Ordering Information

Getting started with the Poppa is as easy as 1, 2, 3

1. Outfit all monitors in the facility with the Poppa.
2. Purchase cuffs with the Poppa connector installed.
3. Easily connect to any monitor.

Connector	Pack Size	Product Number
Poppa 1-Tube Monitor Connector	10/pk	BPC-P001SF
	100/pk	BPC-P001SFB
Poppa 2-Tube Monitor Connector	10/pk	BPC-P002DF
	100/pk	BPC-P002DFB



Poppa Connectors Pre-installed on NIBP Cuffs

All cuffs with Poppa connectors are configured as 1 Tube cuffs and come with the Poppa cuff side connector pre-installed.

CUFF STYLE	Cuff Size/Range		Color Code	Product Number	Cuffs per box
ComfortCheck™ Soft Fabric, Disposable BP Cuffs	Child	11 - 22 cm	GREEN	ZCC4401P	10/bx
	Small Adult	16 - 28 cm	LIGHT BLUE	ZCC4301P	10/bx
	Adult	26 - 33 cm	DARK BLUE	ZCC4201P	10/bx
	Adult Long	26 - 33 cm	DARK BLUE	ZCC4221P	10/bx
	Large Adult	33 - 41 cm	RED	ZCC4101P	10/bx
	Large Adult Long	33 - 41 cm	RED	ZCC4121P	10/bx
	Thigh	39 - 55 cm	BROWN	ZCC4001P	10/bx
StatCheck™ Vinyl, Disposable BP Cuffs	Child	11 - 22 cm	GREEN	ZSC4401P	10/bx
	Small Adult	16 - 28 cm	LIGHT BLUE	ZSC4301P	10/bx
	Adult	26 - 33 cm	DARK BLUE	ZSC4201P	10/bx
	Adult Long	26 - 33 cm	DARK BLUE	ZSC4221P	10/bx
	Large Adult	33 - 41 cm	RED	ZSC4101P	10/bx
	Large Adult Long	33 - 41 cm	RED	ZSC4121P	10/bx
	Thigh	39 - 55 cm	BROWN	ZSC4001P	10/bx
DuraCheck™ Reusable, Bladderless BP Cuffs	Infant	8 - 15 cm	ORANGE	ZDC4901P	5/bx
	Child	11 - 22 cm	GREEN	ZDC4401P	5/bx
	Small Adult	16 - 28 cm	LIGHT BLUE	ZDC4301P	5/bx
	Adult	26 - 33 cm	DARK BLUE	ZDC4201P	5/bx
	Adult Long	26 - 33 cm	DARK BLUE	ZDC4221P	5/bx
	Large Adult	33 - 41 cm	RED	ZDC4101P	5/bx
	Large Adult Long	33 - 41 cm	RED	ZDC4121P	5/bx
Thigh	39 - 55 cm	DARK BLUE	ZDC4001P	5/bx	
VitalCheck™ Reusable Bladder BP Cuffs	Infant	9 - 15 cm	BROWN	ZVC4901P	5/bx
	Child	14 - 21.5 cm	GREEN	ZVC4401P	5/bx
	Small Adult	20.5 - 28.5 cm	LIGHT BLUE	ZVC4301P	5/bx
	Adult	27 - 35 cm	NAVY BLUE	ZVC4201P	5/bx
	Adult Long	27 - 35 cm	NAVY BLUE	ZVC4221P	5/bx
	Large Adult	34 - 44 cm	PLUM	ZVC4101P	5/bx
	Large Adult Long	34 - 44 cm	PLUM	ZVC4121P	5/bx
Thigh	42 - 52 cm	BROWN	ZVC4001P	5/bx	

Distributed By: