

# VitalCheck™



## Reusable Blood Pressure Cuffs

For over ten years Zefon International has been designing and manufacturing superior blood pressure cuffs for leading blood pressure measurement instrument manufacturers. With a complete range of sizes, materials and connectors, Zefon brand cuffs represent the best value to bring consistent performance, patient comfort and economy to your medical setting.

## Features

- Reusable; Durable construction for long life and multiple patient use.
- Cost Effective; Bladder and sleeve design allows sleeves to be replaced as they become worn or contaminated.
- Made from Urethane material that is soft and comfortable for the patient.
- Latex Free.
- Meets AAMI and AHA Cuff Sizing Recommendations.
- Color coded cuff sleeve makes it easy to identify cuff size.
- Designed as a direct replacement for OEM cuffs provided by monitor manufacturers.
- Full range of connectors and adapters allow standardization of cuffs across the whole hospital.
- Easy connection to electronic NIBP monitors or manual blood pressure measurement systems.




5350 SW 1st Lane, Ocala, FL 34474 USA  
Tel: +1 3528548080; Fax: +1 3528547480  
www.zefon.com

MM00005 Rev 7 2010-04

EC REP Blue Box Medical Limited  
Unit 29, New Forest Enterprise Centre,  
Rushington Business Park; Chapel Lane  
Totton, Southampton, SO40 9LA United Kingdom  
Phone: +023 8066 9000; Fax: +023 8066 9002

# Ordering Information

## Single Tube Cuffs




CAS 740/750 Criticare 508 Datascope Accutor Plus & Passport 2 Philips/HP Siemens Spacelabs	BCI, Colin, Datascope Passport XG & Expert Nihon Kodon Lifescopy 9 & 12 Welch Allyn Vital Signs/Atlas (Before July 2000)	Welch Allyn Vital Signs/Atlas (after July 2000)	Welch Allyn Protocol Propaq

Cuff Size	Color Code	Limb Circumference	1-Tube Bayonet Connector	1-Tube Female Luer Connector	1-Tube Screw Type Connector	1-Tube Female Twist Lock Connector
Infant	Brown	9-15cm	ZVC4901HP	ZVC4901F	ZVC4901D	ZVC4901MQF
Child	Green	14-21.5cm	ZVCC4401HP	ZVCC4401F	ZVCC4401D	ZVCC4401MQF
Small Adult	Light Blue	20.5-28.5cm	ZVCC4301HP	ZVCC4301F	ZVCC4301D	ZVCC4301MQF
Adult	Navy Blue	27-35cm	ZVCC4201HP	ZVCC4201F	ZVCC4201D	ZVCC4201MQF
Adult Long	Navy Blue	27-35cm	ZVCC4221HP	ZVCC4221F	ZVCC4221D	ZVCC4221MQF
Large Adult	Plum	34-44cm	ZVCC4101HP	ZVCC4101F	ZVCC4101D	ZVCC4101MQF
Large Adult Long	Plum	34-44cm	ZVCC4121HP	ZVCC4121F	ZVCC4121D	ZVCC4121MQF
Thigh	Brown	42-52cm	ZVC4001HP	ZVC4001F	ZVC4001D	ZVC4001MQF

All cuffs are packaged with 5 cuffs per box

## Double Tube Cuffs



Critikon® Dinamap® Criticare 507 Invivo Research, MDE Nihon Kodon Lifescopy L Spacelabs	Datex-Ohmeda	GE Medical Eagle 4000 Tram 100/200/300	Manual

Cuff Size	Color Code	Limb Circumference	2-Tube Screw Type Connector	2-Tube Male/Female Twist Lock Connector	2-Tube Female Twist Lock Connector	2-Tube Manual Inflation with Male Luer
Infant	Brown	9-15cm	ZVC4902D	ZVC4902MQMF	ZVC4902MQF	ZVC4912L
Child	Green	14-21.5cm	ZVCC4402D	ZVCC4402MQMF	ZVCC4402MQF	ZVCC4412L
Small Adult	Light Blue	20.5-28.5cm	ZVCC4302D	ZVCC4302MQMF	ZVCC4302MQF	ZVCC4312L
Adult	Navy Blue	27-35cm	ZVCC4202D	ZVCC4202MQMF	ZVCC4202MQF	ZVCC4212L
Adult Long	Navy Blue	27-35cm	ZVCC4222D	ZVCC4222MQMF	ZVCC4222MQF	ZVCC4232L
Large Adult	Plum	34-44cm	ZVCC4102D	ZVCC4102MQMF	ZVCC4102MQF	ZVCC4112L
Large Adult Long	Plum	34-44cm	ZVCC4122D	ZVCC4122MQMF	ZVCC4122MQF	ZVCC4132L
Thigh	Brown	42-52cm	ZVC4002D	ZVC4002MQMF	ZVC4002MQF	ZVC4012L

All cuffs are packaged with 5 cuffs per box

Distributed By:



5350 SW 1st Lane,  
Ocala, FL 34474 USA  
Tel: +1 3528548080  
www.zefon.com