

2019-2020 PRODUCT CATALOG



International Edition



Zefon

INTERNATIONAL

The Professionals Choice for Air Sampling Equipment

LABORATORY & SAFETY SUPPLIES



Acetone Vaporizers and Clearing Supplies: 84-85
Microscope Slides & Cover Glass: 86-87
Sample Collection and Storage Bags: 88
Respirator Fit Test Kits: 89-91

Small Wonder Acetone Vaporizer

The Wonder Makers Small Wonder Acetone Vaporizer was developed as an economical and quick way to clear asbestos PCM slides. Specifically designed to meet NIOSH 7400, 7402, or ORM requirements for preparation of PCM samples.



Features:

- Heat control and thermometer built into vaporizer unit
- Padded carrying case protects equipment and supplies
- Detailed instruction manual
- Enough supplies for the analysis of 1/2 gross of slides
- Lightweight but solidly built
- Designed in strict compliance with the NIOSH 7400, 7402, or ORM requirements for preparation of microscope slides for PCM analysis

Specifications

DIMENSIONS	4-7/8" L x 4-1/2" W x 2-3/4" H (Small Wonder vaporizer) 13" L x 11" W x 4" H (Complete Analysis Kit)
MATERIALS	Plastic case, ABS plastic cover on vaporizer, clear plastic bottles, metal and plastic tools.
WEIGHT	2-1/2 lbs. (Small Wonder vaporizer) 8 lbs. (Complete Analysis Kit)
FINISH	Black vaporizer cover with gray cord. Case color may vary.



Complete Kit Shown
Includes everything needed to clear slides for PCM analysis

PRODUCT	DESCRIPTION
WM-42	Small Wonder Analysis Kit Includes acetone vaporizer, 1 box microscope slides, 1/2oz Triacetin, 2oz Acetone, 1 microdropper, 1 syringe, 1 scalpel, tweezers, padded carrying case and instruction manual.
WM-44	Acetone reagent, 2oz, 99% pure
WM-45	Triacetin reagent, 1/2oz, 99% pure
WM-47	1cc Syringe
WM-48	Rounded Edge Scalpel
WM-49	Microdropper



Perm-O-Fix Acetone Vaporizer

A simple, easy to use, acetone vaporizer for use in preparing asbestos PCM samples as per NIOSH 7400.

PRODUCT	DESCRIPTION
ZA0190	Perm-O-Fix Acetone Vaporizer, 120V
ZA0191	Perm-O-Fix Acetone Vaporizer, 220V



HSL/ULO (HSE/NPL) Phase Shift Test Slide

PRODUCT	DESCRIPTION
ZA0210	HSE/NPL Phase Shift Test Slide, Green Band 5
ZA0215	HSE/NPL Phase Shift Test Slide, Red Band 4

Clearing Solution Supplies

PRODUCT	DESCRIPTION
ZA0159	MCE Media for Clearing Solution, 0.5g sheets, 5/pk
ZA0167	Acetone Reagent, ACS Grade, 30ml Bottle
ZA0162	Triacetin, 30 ml Bottle
ZA0164	Dropper Bottle, 1 oz
ZA0165	Amber Bottle, 1 oz
ZA0166	Amber Bottle, 2 oz



Microscope Slides



Premiere® Brand 9100 Series Microscope Slides

Premiere® Brand 9100 Series Microscope Slides are made of corrosion resistant glass with smoothly polished edges. Precleaned and suitable for a wide range of applications. Sized 3" x 1" x approximately 1mm thick.

PRODUCT	DESCRIPTION	QUANTITY
ZA0260	Plain Clear 3" x 1" Microscope Slides	1 Gross (144 slides) 10 Gross (1,440 slides)
ZA0261	Frosted 3" x 1" Microscope Slides	1 Gross (144 slides) 10 Gross (1,440 slides)



Grafc® Brand Microscope Slides

Grafc® Brand Microscope Slides are exceptionally clear, water-white, extra hard glass. Rounded corners, ground and polished edges. Packaged in boxes with paper between each slide.

PRODUCT	DESCRIPTION	QUANTITY
ZA0270	Plain Clear 3" x 1" Microscope Slides	1 Gross (144 slides) 10 Gross (1,440 slides)
ZA0271	Frosted 3" x 1" Microscope Slides	1 Gross (144 slides) 10 Gross (1,440 slides)

Microscope Slide Storage Boxes & Mailers



PRODUCT	DESCRIPTION
CA-MSL-07	Single Slide Mailer, Plastic
CA-MSL-08	Double Slide Mailer, Plastic
CA-MSL-15	5 Slide Holder, Plastic
CA-97-0025	25 Slide Storage Box, White
CA-97-0055	50 Slide Storage Box, White
CA-97-0125	100 Slide Storage Box, White

Microscope Cover Glass

- Packaged in 1 oz boxes that include approximately 100 slides.
- #1.0 Thickness cover glass measures approximately .15mm thick
- #1.5 Thickness cover glass measures approximately .17mm thick

Platinum Line® Cover Glass

PRODUCT	DESCRIPTION
ZA0280	Cover Glass 18mm x 18mm - #1.0 Thickness, 1 oz bx
ZA0281	Cover Glass 18mm x 18mm - #1.5 Thickness, 1 oz bx
ZA0282	Cover Glass 22mm x 22mm - #1.0 Thickness, 1 oz bx
ZA0283	Cover Glass 22mm x 22mm - #1.5 Thickness, 1 oz bx
ZA0284	Cover Glass 22mm x 40mm - #1.0 Thickness, 1 oz bx



Chase Brand Cover Glass

PRODUCT	DESCRIPTION
ZA0290	Cover Glass 22mm x 22mm - #1.0 Thickness, 1 oz bx
ZA0291	Cover Glass 22mm x 22mm - #1.5 Thickness, 1 oz bx
ZA0292	Cover Glass 22mm x 30mm - #1.0 Thickness, 1 oz bx
ZA0293	Cover Glass 22mm x 30mm - #1.5 Thickness, 1 oz bx
ZA0294	Cover Glass 22mm x 40mm - #1.0 Thickness, 1 oz bx
ZA0295	Cover Glass 22mm x 40mm - #1.5 Thickness, 1 oz bx



Cassette Opener

Accessories

PRODUCT	DESCRIPTION
ZA0046	Cassette Opener for 25mm, 37mm and Air-O-Cell® Cassettes
ZA0049	Cassette Opener for 47mm Cassettes
CA-9410	Scalpel, Disposable Plastic with #10 Stainless Blade, 10/bx
CA-9303	Scalpel Handle
CA-9310	Carbon Steel Blades, #10, 100/bx
ZA0053	Kimwipes®, 4.5" x 8.5", 280/bx
ZA0052	Alcohol prep pads, 200/bx
ZA0056	Via-Cell® Linking Solution, Sterile, 15ml
ZA0058	Via-Cell® Blue Stain, 15ml (For staining the sample for microscopic analysis when culture is not required)

Whirl-Pak® Bags



Nasco Whirl-Pak® sample collection bags provide an easy to use leakproof closure for collection of a wide variety of sample materials.

Features:

- Models with write-on strip make sample identification easier and more convenient by working with a wide variety of writing instruments including ballpoint pens, pencil, or felt-tip markers.
- Puncture proof tabs • Leakproof Closure • Guaranteed Sterile (sterilized after manufacturing)

Whirl-Pak® Bags Without Write-On Area

PRODUCT	DESCRIPTION	U/M
B01009WA	3"w x 5"l, 2.25mil, 2oz capacity	500/bx
B00679WA	3"w x 7.25"l, 2.25mil, 4oz capacity	500/bx
B01018WA	5.25"w x 7.5"l, 3.0mil, 13oz capacity	500/bx
B00736WA	4.5"w x 9"l, 2.5mil, 18oz capacity	500/bx
B01020WA	6"w x 9"l, 3.0mil, 24oz capacity	500/bx
B01027WA	6"w x 15"l, 3.0mil, 42oz capacity	500/bx



The white, write-on strip makes sample identification easier and more convenient. It's easy to write on this white strip using an ordinary ballpoint pen, pencil, or felt-tip marker. The paper-like finish virtually eliminates ink smearing, and provides a surface that will accept a variety of writing instruments.

Whirl-Pak® Bags With Write-On Area

PRODUCT	DESCRIPTION	U/M
B01064WA	3"w x 5"l, 2.25mil, 2oz capacity, includes write on area	500/bx
B01062WA	3"w x 7.25"l, 2.25mil, 4oz capacity, includes write on area	500/bx
B01489WA	3.5"w x 7"l, 2.5mil, 7oz capacity, includes write on area	500/bx
B01065WA	4.25"w x 9"l, 2.5mil, 18oz capacity, includes write on area	500/bx
B01297WA	6"w x 9"l, 2.5mil, 24oz capacity, includes write on area	500/bx
B01195WA	7.5"w x 12"l, 4.0mil, 55oz capacity, includes write on area	500/bx

Zip-Style Bags

These zip style reclosable bags with white write-on area are perfect for labeling custom orders, indicating manufacturing dates, assembly instructions or warnings, as well as labeling medium or heavy weight parts and samples. White-Block area accepts marker, rubber stamp, grease pencil, or pen. Side-seal construction strengthens the reclosable top. Meets FDA/USDA specifications.



PRODUCT	DESCRIPTION	U/M
ZB1	2.5" x 3", 2mil with write on area	100/pk
ZB5	3" x 5", 2mil with write on area	100/pk
ZB4	3" x 5", 4mil with write on area	100/pk
ZB6	9" x 12", 2mil with write on area	100/pk
ZB2	9" x 12", 4mil with write on area	100/pk

Saccharin Respirator Fit Test Kit

The Saccharin Fit Test Kit offers OSHA-compliant testing of reusable and disposable dust and mist respirators. When the sweet taste of saccharin is detected, the respirator has failed its fit test. It includes polycoated hood, two nebulizers, sensitivity and fit test solutions, and instructions.



PRODUCT	DESCRIPTION	U/M
AL-2040	Saccharin Respirator Fit Test Kit	Each
AL-2040-11K	Saccharin Sensitivity Solution Tubes	6/bx
AL-2040-12K	Saccharin Test Solution Tubes	6/bx
AL-2040-04	Replacement Fit Testing Hood	Each
AL-2040-08N	Nebulizer with Rubber Bulb	Each

Banana Oil Respirator Fit Test Kit

Iso Amyl Acetate (banana oil) offers the oldest method of testing respirator fit. Its strong, readily identifiable odor will be recognized if respirators don't obtain a fit. For quick fit testing the standard kit contains two 1 oz. bottles of test agent, vinyl test enclosure, measuring pipette, and paper towels in a handy carry bag.



PRODUCT	DESCRIPTION	U/M
AL-0203	Banana Oil Respirator Fit Test Kit	Each
AL-0201	Banana Replacement Ampules	10/bx
AL-0203-01	Banana Oil Solution	Each

VeriFit® Irritant Smoke Generator for Respiratory Fit Testing

Patented all-in-one design is convenient and complete; no components to lose. Its small integrated bellows pump provides the optimum amount of smoke for fit testing yet minimizes the risk of overexposure. The plastic smoke generator eliminates danger from exposed broken glass for safer operation and disposal.

Features:

- Integrated pump provides correct amount of smoke each time
- All-in-one design - No parts to lose or replace
- Plastic Design - No glass to break



NEXTTEQ
BETTER TECHNOLOGY. BETTER QUALITY.
BETTER VALUE.

Respirator fit testing is accomplished in three simple steps:

1. Bend the smoke generator to break the enclosed ampoule.
2. Remove the cap
3. Wait 5 seconds and compress bellows

PRODUCT	DESCRIPTION	U/M
90095	VeriFit® Irritant Smoke Generator Kit	6/pk
50811000-310N	VeriFit® Irritant Smoke Generator Kit	10/pk

Bitrex® Respirator Fit Test Kit

The Bitrex® Fit Test Kit offers qualitative OSHA-compliant fit testing for disposable and reusable dust and mist respirators. Bitrex creates an unmistakable bitter taste that indicates break through in the respirator fit. Bitrex has been found safer than saccharin tests while producing identical results. Disposable glass ampoules make pouring solution into nebulizers easier and reduces the chance of contamination.

Features:

- For testing the facepiece-to-face seal
- Fast, easy method for performing qualitative fit testing
- Use with any particulate respirator or any gas/vapor respirator with particulate prefilter
- Helps satisfy respiratory program requirements
- Can be used during respirator training programs
- Helps meet administrative respiratory protection training program regulations as specified in OSHA 29 CFR 1910.134



PRODUCT	DESCRIPTION	U/M
AL-2041	Bitrex® Respirator Fit Test Kit Includes hood, 2 nebulizers, fit test solution, sensitivity solution and instructions	Each
AL-2041-11K	Sensitivity Solution Tubes	6/bx
AL-2041-12K	Test Solution Tubes	6/bx

Smoke Tubes Respirator Fit Test Kit

Field proven as an OSHA accepted method since 1982, irritant smoke qualitative respirator fit testing is an economical, fast, safe and reliable fit testing method. Smoke tubes are used with a wide variety of respirators by firefighters, HAZMAT teams, first responders, the military and many industrial applications.



Kit Includes:

- 6 - Irritant Smoke Tubes
- 1 - Precision Aspirator Bulb Assembly
- 1 - Record Keeping Form
- 1 - Carrying Pouch

Irritant smoke testing is fast, simple, easy-to-perform and economical. You save money because there is no need for nebulizers, electric pumps, mixing jars, batteries, special test masks, mask probes or testing hoods. Each kit has everything you need for over 50 tests without the cost of extra equipment. Plus, you save time because each smoke tube and aspiration bulb pump is pre-calibrated and ready for immediate use. No warm up time or calibration required.

PRODUCT	DESCRIPTION	U/M
AL-500	Smoke Tube Kit Contains 6 smoke tubes, aspirator bulb and carrying case	Each
AL-501	Replacement Smoke Tubes 6/bx	6/bx

Trademarks

All trademarks appearing in this catalog are the property of their respective owners. Below is a list of trademarks used within this catalog and their respective owners.

Zefon: Air-O-Cell, Bio-Pump, DigiCal, EconoGrab, Escort, Escort ELF, Gemini, Gravi-Sert, MoldSnap, Solu-Sert, Via-Cell, Z-Lite, ZePore and **Zefon** are trademarks or registered trademarks of **Zefon International, Inc.**

3M Company: 3M, ChromateCheck and LeadCheck are trademarks or registered trademarks of 3M Company

AP Buck: Allin, Bio-Culture, Elite, and Libra are trademarks of AP Buck Inc.

Buchem: Tenax is a registered trademark of Buchem BV Corporation

C & A: Premiere is a registered trademark of C & A Scientific Co. Inc.

Dupont: Tedlar is a registered trademark of E. I. Du Pont De Nemours and Company

Engineering Solutions: Omniguard is a trademark of Engineering Solutions, Inc.

Environmental Monitoring Systems: Micro5 is a trademark of Environmental Monitoring Systems

GE Healthcare: 934AH is a registered trademark of GE Healthcare/Whatman, Inc.

GF Health: Grafo is a registered trademark of GF Health Products Inc.

Imerys: Chromosorb is a registered trademark of Imerys Minerals California

Johnson Matthey: Bitrex is a registered trademark of Johnson Matthey Public Limited Company

Kimberly Clark: Kimwipes is a registered trademark of Kimberly-Clark Worldwide, Inc.

Mercedes Medical: Platinum Line is a registered trademark of Mercedes Medical Inc.

Microsoft: Excel and Microsoft are registered trademarks of Microsoft Corporation


Nasco: Whirl-Pak is a registered trademark of Nasco/The Aristotle Corporation

Nextteq: Nextteq, VeriAir Flex and VeriFit are trademarks of Nextteq LLC

Pall: Teflo, Tissuquartz, and Zefluor are trademarks of Pall Corporation

Pelican Products: Pelican is a trademark of Pelican Products Inc.

St. Gobain: Tygon is a registered trademark of Saint Gobain Performance Plastic Corporation



**It shouldn't be
such a struggle to
get the products
you need.**

Zefon has the air sampling
products you're looking for to
make your buying experience
easy and hassle free.

Zefon

INTERNATIONAL


The Professionals Choice for Air Sampling Equipment

Contact your local representative for more information | Follow us on



Zefon

INTERNATIONAL



Follow us to stay up-to-date on new products, specials, and company news!



www.zefon.com



Zefon

INTERNATIONAL

The Professionals Choice for Air Sampling Equipment